

There's a life behind every test

The **BD Veritor™ Plus System** enhances point-of-care testing for **SARS-CoV-2**,* Flu A+B, RSV,† and Group A Strep assays

Discover real progress in point-of-care (POC) testing with the **BD Veritor™ Plus System**. Combine reliable results with ease of implementation and simplified workflow options.



*EUA authorized by FDA
†respiratory syncytial virus



BD Veritor™ Plus System

Taking the point-of-care
experience to the **next level**

*The intended use of the BD Veritor™ Plus System for Rapid Detection of SARS-CoV-2 assay only includes those who are suspected of COVID-19 by their health care provider within the first five days of the onset of symptoms.

- This test has not been FDA cleared or approved
- This test has been authorized by FDA under an Emergency Use Authorization (EUA) for use by authorized laboratories
- This test has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens
- This test is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Act, 21 U.S.C. § 360bbb-3(b)(1), unless the authorization is terminated or revoked sooner

At BD, we understand that point-of-care testing presents challenges:

- Complex **manual workflows**
- Lengthy time to results and **delays in diagnosis**
- **Manual documentation**, potentially leading to transcription errors
- Inadequate **training programs** or **reference materials**
- Difficulties in **checking training compliance** or **assessing staff competency**

The BD Veritor™ Plus Analyzer is a portable, easy-to-use instrument that allows you to overcome challenges to providing streamlined POC testing



Offers flexible and adaptable workflow options, allowing you to maximize staff efficiency



Minimizes manual transcription by allowing data storage/download via USB port



Delivers timely results — 15 minutes for SARS-CoV-2,¹ less than 11 minutes for Flu A+B² and RSV,³ and less than 6 minutes for Group A Strep⁴ — minimizing patient wait times



Includes the interactive BD eLearning online training platform, providing a comprehensive curriculum

- **Automatically updated with the latest** resources and references
- **Allows you to review training compliance and manage certificates** of program completion



Addressing the many challenges of point-of-care testing

The **BD Veritor™ Plus System**, together with the BD eLearning online training platform, provide you with:

- Optimal training
- Easy implementation and use
- Workflow flexibility
- Reliable, timely results
- Streamlined sample processing

✓ Optimal training

The BD eLearning online training platform provides a comprehensive experience with:

- Interactive, step-by-step training modules
- Dashboards and reports for easy management of training compliance and completion

✓ Easy implementation and use

The **BD Veritor™ Plus System** simplifies operation by offering:

- Easy-to-use tests, allowing novice and non-lab personnel to confidently help with diagnosis of SARS-CoV-2, Flu A+B, RSV, and Group A Strep
- Intuitive, single-button functionality, reducing training time
- A digital display for easy-to-read results

✓ Workflow flexibility

Select the mode that fits your workflow

- To keep up with high test demand by batching, use the **Analyze Now** mode
 - » Incubate the specimen for the time required for each assay, then insert the test device into the Analyzer; the result is displayed within seconds
- To accomplish other tasks while running tests, use the **Walk Away** mode
 - » Add the specimen to the test device and insert into the Analyzer for automatic incubation and analysis



✓ Reliable, timely results

Proven performance vs reference methods for SARS-CoV-2, Flu A+B, RSV, and Group A Strep

	CLIA Status	PPA [*]	NPA [†]	Time to Result
SARS-CoV-2¹ (vs RT-PCR)	-	84% (95% CI, 67%-93%)	100% (95% CI, 98%-100%)	15

SARS-CoV-2 result in 15 minutes

	CLIA Status	PPA [*]	NPA [†]	Time to Result
Influenza A² (vs PCR)	CLIA-waived	83.6% (76.1%, 89.1%)	97.5% (95.7%, 98.5%)	<11
Influenza B² (vs PCR)	CLIA-waived	81.3% (71.1%, 88.5%)	98.2% (95.7%, 99.3%)	

Flu A+B result in <11 minutes


	CLIA Status	Sensitivity	Specificity	Time to Result
RSV³ (vs viral culture)	CLIA-waived	91.8% (95% CI, 85.9%-95.4%)	93.3% (95% CI, 90.4%-94.4%)	<11
Group A Strep⁴ (vs culture)	CLIA-waived	95.4% (95% CI, 90.3%-97.9%)	95.7% (95% CI, 93.7%-97.1%)	<6

RSV result in <11 minutes; Group A Strep result in <6 minutes

✓ Streamlined sample processing

Easy processing uses prefilled, unitized tubes color coded by reagent or assay type, reducing potential errors.

To process the sample:



1. Insert the sample swab into the tube (after adding Reagent 1 for Group A Strep)
2. Mix the sample with the reagent
3. Dispense sample into the test device
4. Insert test device into the Analyzer:
 - For **Walk Away** mode, insert immediately
 - For **Analyze Now** mode, insert after processing the sample


*positive percent agreement
†negative percent agreement

Personalized solutions for training and quality control

The BD eLearning online training platform offers personalized instruction, allowing staff to quickly become proficient with the **BD Veritor™ Plus System** for smoother implementation

Comprehensive training modules for all available assays (SARS-CoV-2, Flu A+B, RSV, and Group A Strep)

- Includes introduction to the **BD Veritor™ Plus System**
- Each **5-minute module** follows the same structure:



Delivers a personalized, dynamic learning experience

- Allows the assignment of training plans by role, product, or learning type
- Facilitates learning with interactive videos featuring step-by-step instructions and concrete examples, available on any device, any place, any time

Helps track training compliance

- Helps ensure training compliance with easy-to-use reporting and dashboard tools
- Allows convenient management of certificates of program completion

Provides a complete library of resources and references

Updated with the latest FDA-cleared versions of:



The BD Veritor™ InfoScan module facilitates **reliable quality control**

Scanning capabilities facilitate traceability and may reduce hands-on time, transcription errors, and paperwork

Automated documentation functionality includes:

- Operator ID
- Patient/specimen ID
- Kit lot information

Barcode scanner can read the full 128 ASCII character set, including all US keyboard characters

- Data can be stored locally or downloaded via USB cable for improved traceability (USB cable included)
- Results can be printed



Elevate the point-of-care experience for patients and staff

- Easy-to-read display featuring single-button functionality and minimal hands-on time to meet your staff's demands
- Flexible workflow options to adapt to your workplace needs
- Timely, reliable results (15 minutes for SARS-CoV-2¹; <11 minutes for Flu A+B² and RSV³; <6 minutes for Group A Strep⁴), minimizing wait times
- BD eLearning online training platform that simplifies staff instruction and system implementation
- Improved quality control and traceability enabled by the optional BD Veritor™ InfoScan module

Ordering information	Cat no.	Qty
BD Veritor™ Plus System SARS-CoV-2	256082	30 tests
BD Veritor™ Plus System Flu A+B CLIA-waived kit	256045	30 tests
BD Veritor™ Plus System Flu A+B (moderately complex)	256041	30 tests
BD Veritor™ Plus System RSV CLIA-waived kit	256038	30 tests
BD Veritor™ Plus System RSV (moderately complex)	256042	30 tests
BD Veritor™ Plus System Group A Strep CLIA-waived kit	256040	30 tests
BD Veritor™ Plus Analyzer	256066	1 unit
BD Veritor™ InfoScan module	256068	1 unit
USB printer cable	443907	1 unit

References: 1. BD Veritor™ Plus System for Rapid Detection of SARS-CoV-2. Package insert. Becton, Dickinson and Company. 2. BD Veritor™ Plus System for Rapid Detection of Flu A+B. Package insert. Becton, Dickinson and Company. 3. BD Veritor™ Plus System for Rapid Detection of Respiratory Syncytial Virus (RSV). Package insert. Becton, Dickinson and Company. 4. BD Veritor™ Plus System for Rapid Detection of Group A Strep. Package insert. Becton, Dickinson and Company.

BD Life Sciences, 7 Loveton Circle, Sparks, MD 21152-0999 USA
800.638.8663

BDVeritor.com

BD, the BD Logo and Veritor are trademarks of Becton, Dickinson and Company or its affiliates.
© 2020 BD. All rights reserved. 670-US-1120 December 2020

